



Genesys Spine Apache™ TLIF System



EMERGO EUROPE

Molenstraat 15

2513-BH, The Hague


The Netherlands



1250 Capital of Texas Highway S.
Building 3, Suite 600
Austin, TX 78746

LC-009 REV B

The Genesys Spine Apache™ TLIF System consists of precision instruments and PEEK Optima® implants for aiding in fusion of the lumbar spine. The combination of PEEK Optima® and tantalum markers allows surgeons to have an enhanced view of implant placement during surgery and vertebral fusion post surgery.



HEIGHTS*

(*1mm increments)

6mm - 15mm

SPECIFICATIONS

Length = 35mm

Width (Anterior arc) = 14mm

Available in 8° Lordotic or Parallel

Genesys
spine™

*FDA approved for use as an IBFD and VBR

1250 Capital of Texas Highway S.
Building 3, Suite 600
Austin, TX 78746