



---

## Genesys Spine Apache™ PLIF System

---



EMERGO EUROPE

Molenstraat 15

2513-BH, The Hague


The Netherlands

Genesys  
spine™

1250 Capital of Texas Highway S.  
Building 3, Suite 600  
Austin, TX 78746

LC-010 REV B

The Genesys Spine Apache™ PLIF System consists of precision instruments and PEEK Optima® implants for aiding in fusion of the lumbar spine. The combination of PEEK Optima® and tantalum markers allows surgeons to have an enhanced view of implant placement during surgery and vertebral fusion post surgery.



## HEIGHTS\*

(\*1mm increments)

6mm - 15mm

## SPECIFICATIONS

Available in 25mm and 30mm lengths

9.35mm width

Available in Lordotic, Parallel, and Convex

**Genesys**  
spine™

1250 Capital of Texas Highway S.  
Building 3, Suite 600  
Austin, TX 78746

\*FDA approved for use as an IBFD and VBR