



---

## Genesys Spine Apache™ Cervical Interbody System


---

Genesys  
spine™

1250 Capital of Texas Highway S.  
Building 3, Suite 600  
Austin, TX 78746

LC-011 REV B

The Genesys Spine Apache™ Cervical Interbody System consists of precision instruments and PEEK Optima® implants for aiding in fusion of the cervical spine. The combination of PEEK Optima® and tantalum markers allows surgeons to have an enhanced view of implant placement during surgery and vertebral fusion post surgery.



## HEIGHTS\*

(\*1mm increments)

5mm - 12mm

## SPECIFICATIONS

Small: 12mm x 14mm

Large: 15mm x 17mm

Available in 6° Lordotic, Parallel, and Convex



EMERGO EUROPE

Molenstraat 15

2513-BH, The Hague

The Netherlands

**Genesys**  
spine™

1250 Capital of Texas Highway S.  
Building 3, Suite 600  
Austin, TX 78746