

Genesys Spine Apache™ ALIF System



EMERGO EUROPE

Molenstraat 15
2513-BH, The Hague
The Netherlands



1250 Capital of Texas Highway S.
Building 3, Suite 600
Austin, TX 78746

LC-012 REV B


Genesys Spine Apache™ ALIF System

p: +1.512.381.7070

f: +1.512.381.7076

info@genesyspine.com

The Genesys Spine Apache™ ALIF System consists of precision instruments and PEEK Optima® implants for aiding in fusion of the lumbar spine. The combination of PEEK Optima® and tantalum markers allows surgeons to have an enhanced view of implant placement during surgery and vertebral fusion post surgery.



HEIGHTS*

(*2mm increments)

8mm - 20mm
(6° Lordotic)

10mm - 20mm
(12° Lordotic)

SPECIFICATIONS

Length = 25mm

Width = 35mm

Available in 6° and 12° Lordotic

Genesys
spine™

*FDA approved for use as an IBFD and VBR

1250 Capital of Texas Highway S.
Building 3, Suite 600
Austin, TX 78746

Refer to the product insert for detailed indications/contraindications, warnings/precautions and possible adverse effects. To obtain labeling limitations, surgical technique manuals, and/or more information regarding Genesys products, contact your local sales representative or Genesys customer service, at +1 512 381-7070.