


Apache® Cervical Interbody System

Apache® Cervical Interbody System

The Apache® Cervical Interbody System consists of precision instruments and PEEK Optima® implants for aiding in fusion of the cervical spine. The combination of Peek OPTIMA® and tantalum markers allows surgeons to have an enhanced view of implant placement during surgery and vertebral fusion post surgery.



SPECIFICATIONS

Small: 12mm x 14mm
15mm x 19mm

Medium: 14mm x 17mm

Large: 15mm x 17mm
12mm x 15mm

HEIGHTS*

(*1mm increments)

5mm - 12mm

OPTIONS

5° and 7° Lordotic, Parallel and Convex



Lordotic



Parallel



Convex



Patent: www.genesyspine.com/patents

LC-011 Rev D